


# L11 DNA Assembly

Instructor: Prof. K. Solomon Ph.D.  
Assistant Professor  
Agricultural & Biological Engineering  
Laboratory for Renewable  
Resources Engineering

Fall 2018



## Recall...



2

## This lecture....

- Vector features
- Assembly basics
- BioBrick™ Standard

3

## Design fundamentally rooted in DNA


How do we deliver our constructs?

1. Plasmids – circular replicating DNA
2. Host genome
  - Integrated in natural chromosome
  - Synthetic chromosomes

DNA must be replicated for it to be useful

4

## Plasmids contain at least 3 key components



idtdna.com

5

**AB Resistance or Selectable Marker** – provides pressure for the cell to retain plasmid

**Origin (ori)** – site of replication initiation/essential for replication

**MCS** – region where you place your DNA design

Plasmid nomenclature  
pXXX e.g. pUC19, pBAD24

## Selection markers

Plasmids cost resources to maintain

- selection gives the cell reason to propagate your construct

**Bacteria:** Selection is typically with an antibiotic that inhibits or prevents growth

- ampicillin – bacteriostatic, inhibits cell wall synthesis (*bla*)
- chloramphenicol – inhibits protein synthesis/23s rRNA (*cat*)
- kanamycin – inhibits protein synthesis/30S subunit (*neo*)
- tetracycline – inhibits protein synthesis/rRNA:ribosome (*tetA*)

**Eukaryotes (yeast):** Selection is typically with a nutrient auxotrophy

- Uracil – *URA3* – needed for RNA synthesis
- Leucine – *LEU2* – needed for protein synthesis
- Histidine – *HIS3* – needed for protein synthesis

*can also be used for counterselection – e.g. 5-FOA is assimilated by URA3 to make a toxic compound*

6

## Origins of replication

Allows replication and **sets plasmid copy number** in the cells

*How does plasmid copy number affect gene expression?*

### Yeast

- CEN ~ 1 copy per cell (centromeric)
- 2 micron ~ 50 copies per cell

\*Yeast integrating plasmids have no ORI as they must be linearized and integrated into host for replication

Different yeast plasmids may be combined as long as each have a distinct selectable marker

7

## Origins of replication


### Bacteria

Common Vectors	Copy Number*	ORI	Incompatibility Group
pUC	~500-700	pMB1 (derivative)	A
pBR322	~15-20	pMB1	A
pET	~15-20	pBR322	A
pGEX	~15-20	pBR322	A
pColE1	~15-20	ColE1	A
pR6K	~15-20	R6K*	C
pACYC	~10	p15A	B
pSC101	~5	pSC101	C
pBluescript	~300-500	ColE1 (derivative) and F1**	A
pGEM	~300-500	pUC and F1**	A

*Only plasmids from different incompatibility groups with unique selection markers may be combined*

8


### Incompatibility groups share replication mechanisms: ColE1 type A plasmids in *E. coli*



Klumpp S (2011) Growth-Rate Dependence Reveals Design Principles of Plasmid Copy Number Control. PLoS ONE 6(5): e20403.  
doi:10.1371/journal.pone.0020403  
<http://journals.plos.org/plosone/article?id=doi:10.1371/journal.pone.0020403>



### Copy number set by Negative Feedback: ColE1 type plasmids in *E. coli* (Feedback via antisense RNA)



Copy number set by relative expression of RNAII and RNAI and their interactions

e.g. pUC vs pMB1  
Rop<sup>-</sup> vs Rop<sup>+</sup>  
~300 vs ~20

Other plasmid incompatibility groups have other mechanisms

<https://en.wikipedia.org/wiki/ColE1>

## Plasmid types

**Expression plasmids** will typically include:

- promoters:  $P_{T7}$ ,  $P_{Tet}$ ,  $P_{lac}$ ,  $P_{BAD}$ , GALP
- Regulatory proteins: *lacI*, *tetR*, *araC*
- Terminators:  $T7_{term}$ , *rrnB*

**Shuttle vectors** will have features for multiple species to "shuttle" DNA between them (e.g. *ori* and *selecton*)




- E.g. pRS314 – bacterial features for cloning, yeast features for expression

11

## Methods for DNA assembly

How do we put together our circuits?

1) Chemical addition (oligonucleotide/primer synthesis)



12

## Methods for DNA assembly

How do we put together our circuits?

- 2) Homologous Recombination  
e.g. CRISPR, Recombineering (Genome integration)

**“Nucleotide sequences are exchanged between similar sequences of DNA”**




13

## Methods for DNA assembly

How do we put together our circuits?

- 3) Anneal + Polymerization (Gibson CBA, SLIC, commercial DNA synthesis)

*Exploits the ability of single stranded DNA to anneal (stick) to complementary strands*



Wikipedia.com


14

## Methods for DNA assembly

How do we put together our circuits?

- 4) Digestion and Ligation

*Exploits restriction enzymes to create ‘sticky’ overhangs that can be annealed to join pieces of DNA*



idtdna.com


15

## Restriction Enzymes


**Most** restriction enzymes recognize and cut unique 6nt **palindromic** sequences

→ 5' – 3' on either strand is exactly the same

EcoRI



Res often make staggered cuts yielding overhangs



16

## Digestion + Ligation

DNA cut with the same RE have the same overhang and anneal through complementary H-bonds (palindromic!)

H-bond immobilizes fragments long enough for permanent ligation by **ligase**

*How do we standardize this?*

17

## BioBrick™ Standard

Standardizing on a *single* RE doesn't allow you to add assemble parts into devices and other higher order structures

*Why?*

Biobricks use 4 REs:


EcoRI	5'...GAATTC...3'	5'...A <sup>↓</sup> CTAGT...3'	SpeI
	3'...CTTAAG...5'	3'...TGATC <sup>↓</sup> A...5'	
PstI	5'...CTGCAG...3'	5'...T <sup>↓</sup> CTAGA...3'	XbaI
	3'...GACGTC...5'	3'...AGATC <sup>↓</sup> T...5'	

**SpeI and XbaI have complementary overhangs that destroy the restriction site after ligation**

18

## BioBrick™ Assembly

All parts have a prefix (EcoRI and XbaI, EX) and a suffix (SpeI and PstI, SP)



*How would you do this if you wanted to put C0010 in front of B0034?*

Ligation site has a **scar**: ACTAGA  
TGATCT

SpeI XbaI

19

## BioBrick™ Assembly

Advantages	Disadvantages
<ul style="list-style-type: none"> <li>- idempotent/scalable</li> <li>- Easily automatable</li> <li>- Robust</li> <li>- Only 4 REs used</li> </ul>	<ul style="list-style-type: none"> <li>- Need to remove internal Biobrick RE sites from all parts</li> <li>- Have a fixed scar (not good for protein fusions)</li> <li>- Can't add parts between parts</li> </ul>

20

### **Next time....**

---

- More advanced high-throughput assembly techniques

21